

## EXPERIENCE STRENGTH & STABILITY

PRECISION SPINE  
**REFORM<sup>®</sup> Ti**  
TITANIUM PEDICLE SCREW SYSTEM

Low profile titanium tulip  
**REDUCES RISK OF FACET IMPINGEMENT**

Titanium square threads  
**REDUCE RISK OF CROSS-THREADING**

30° of angulation in all planes  
**ALLOWS FOR INTRAOPERATIVE FLEXIBILITY**

Triple lead thread  
**INCREASES SPEED & EFFICIENCY OF SCREW DELIVERY  
AND REDUCES FATIGUE**



Enlarged T25 drive feature  
**PROVIDES STABILITY DURING INSERTION**

Minimally tapered screw thread  
**LOWERS SCREW INSERTION TORQUE  
WHILE OPTIMIZING PULL OUT  
STRENGTH**

Self starting aggressive screw tip  
**OPTIMIZES INITIAL BONE  
ENGAGEMENT**

**PRECISION SPINE<sup>®</sup>**  
*Discover the Difference*



## SYSTEM TRAYS

P/N	Description
39-BK-070X	Reform Ti Standard Implants
39-BK-0721	Reform Ti Standard Instruments Tray 1
39-BK-0202	Reform Ti Standard Instruments Tray 2



## BONE SCREWS

### Ø 4.5mm SCREWS

P/N	Description
39-TP-4525*	4.5 x 25mm - Reform Ti Screw
39-TP-4530*	4.5 x 30mm - Reform Ti Screw
39-TP-4535*	4.5 x 35mm - Reform Ti Screw
39-TP-4540*	4.5 x 40mm - Reform Ti Screw
39-TP-4545*	4.5 x 45mm - Reform Ti Screw

### Ø 5.5mm SCREWS

P/N	Description
39-TP-5530	5.5 x 30mm - Reform Ti Screw
39-TP-5535	5.5 x 35mm - Reform Ti Screw
39-TP-5540	5.5 x 40mm - Reform Ti Screw
39-TP-5545	5.5 x 45mm - Reform Ti Screw
39-TP-5550	5.5 x 50mm - Reform Ti Screw

### Ø 6.5mm SCREWS

P/N	Description
39-TP-6530	6.5 x 30mm - Reform Ti Screw
39-TP-6535	6.5 x 35mm - Reform Ti Screw
39-TP-6540	6.5 x 40mm - Reform Ti Screw
39-TP-6545	6.5 x 45mm - Reform Ti Screw
39-TP-6550	6.5 x 50mm - Reform Ti Screw
39-TP-6555	6.5 x 55mm - Reform Ti Screw

### Ø 7.5mm SCREWS

P/N	Description
39-TP-7535	7.5 x 35mm - Reform Ti Screw
39-TP-7540	7.5 x 40mm - Reform Ti Screw
39-TP-7545	7.5 x 45mm - Reform Ti Screw
39-TP-7550	7.5 x 50mm - Reform Ti Screw
39-TP-7555	7.5 x 55mm - Reform Ti Screw

### Ø 8.5mm SCREWS

P/N	Description
39-TP-8540*	8.5 x 40mm - Reform Ti Screw
39-TP-8545*	8.5 x 45mm - Reform Ti Screw
39-TP-8550*	8.5 x 50mm - Reform Ti Screw
39-TP-8555*	8.5 x 55mm - Reform Ti Screw
39-TP-8560*	8.5 x 60mm - Reform Ti Screw
39-TP-8570*	8.5 x 70mm - Reform Ti Screw
39-TP-8580*	8.5 x 80mm - Reform Ti Screw

### Ø 9.5mm SCREWS

P/N	Description
39-TP-9560*	9.5 x 60mm - Reform Ti Screw
39-TP-9570*	9.5 x 70mm - Reform Ti Screw
39-TP-9580*	9.5 x 80mm - Reform Ti Screw

\* Per Request Only

# PRECISION SPINE REFORM<sup>®</sup>Ti TITANIUM PEDICLE SCREW SYSTEM

## LOCKING CAP, CONNECTORS & RODS

### LOCKING CAP

P/N	Description
39-LS-0100	Cap Screw (Lock-Screw)



### CROSS CONNECTORS

P/N	Description
39-CC-0035	Cross Connector 35mm
39-CC-0040	Cross Connector 40mm
39-CC-0048	Cross Connector 48mm
39-CC-0066	Cross Connector 66mm



### Ø 5.5mm LORDOTIC Ti RODS

P/N	Description
39-LT-5035	Lordotic Ti Rod, x 35mm
39-LT-5040	Lordotic Ti Rod, x 40mm
39-LT-5045	Lordotic Ti Rod, x 45mm
39-LT-5050	Lordotic Ti Rod, x 50mm
39-LT-5055	Lordotic Ti Rod, x 55mm
39-LT-5060	Lordotic Ti Rod, x 60mm
39-LT-5065	Lordotic Ti Rod, x 65mm
39-LT-5070	Lordotic Ti Rod, x 70mm
39-LT-5075	Lordotic Ti Rod, x 75mm
39-LT-5080	Lordotic Ti Rod, x 80mm
39-LT-5090	Lordotic Ti Rod, x 90mm
39-LT-5100	Lordotic Ti Rod, x 100mm
39-LT-5110	Lordotic Ti Rod, x 110mm
39-LT-5120	Lordotic Ti Rod, x 120mm



### Ø 5.5mm STRAIGHT Ti RODS

P/N	Description
39-ST-5080	Straight Ti Rod, x 80mm
39-ST-5100	Straight Ti Rod, x 100mm
39-ST-5120	Straight Ti Rod, x 120mm
39-ST-5200	Straight Ti Rod, x 200mm
39-ST-5400	Straight Ti Rod, x 400mm



PRECISION SPINE<sup>®</sup>  
www.precisionspineinc.com

