


REDEFINING RESULTS

45° max screw angulation
MAXIMIZES BONE PURCHASE & STABILITY

Large graft window
OPTIMIZES BIOLOGICAL
INGROWTH TO
FACILITATE FUSION

One-step locking mechanism
PREVENTS BONE SCREW BACKOUT

PRECISION SPINE
DAKOTA ACDF™
STANDALONE CERVICAL SYSTEM



Dual thread screws
MAXIMIZE CORTICAL
AND CANCELLOUS
BONE INTERFACE

7° Lordosis, 2 Footprints and
6-12mm heights (1mm incr.)
ACCOMMODATE VARYING
PATIENT ANATOMY



Integrated
Plate and Cage
Interbody

PRECISION SPINE®
Discover the Difference



Instruments

Item No.	Description	Qty
71-CH-0001	AO Modular Handles	2
37-IN-0060	Screw/Locking Driver, T8	2
71-TN-7406	Trial, 14mm x 12mm, 5/6mm, 7°	1
71-TN-7408	Trial, 14mm x 12mm, 7/8mm, 7°	1
71-TN-7410	Trial, 14mm x 12mm, 9/10mm, 7°	1
71-TN-7412	Trial, 14mm x 12mm, 11/12mm, 7°	1
71-TN-7606	Trial, 16mm x 14mm, 5/6mm, 7°	1
71-TN-7608	Trial, 16mm x 14mm, 7/8mm, 7°	1
71-TN-7610	Trial, 16mm x 14mm, 9/10mm, 7°	1
71-TN-7612	Trial, 16mm x 14mm, 11/12mm, 7°	1
71-IN-0030	Drill Guide	1
71-SP-0212	AO Straight Drill - 12mm x 1.7mm	1
71-SP-0214	AO Straight Drill - 14mm x 1.7mm	1
71-SP-0215	AO Straight Drill - 15mm x 1.7mm	1
71-IN-0010	Straight Awl, 9mm x 1.7mm	1
71-IN-0720	Interbody Insert	2
37-IN-0050	Tamp	1
37-IN-0084	Impaction/Extraction Mallet	1
Angled Instruments		
71-IN-0020	Angled Awl, 40°, 9mm x 1.7mm	1
71-IN-0410	AO Angled Drill - 12mm x 1.7mm	1
71-IN-0400	AO Angled Driver, 40°, T8	1
71-IN-0407	Counter Torque Handle	1
By Request Instruments		
71-IN-0415	AO Angled Driver, Short Tip, 40°, T8 (NOT COMPATIBLE with Interbody Insert)	1
71-BK-0110	Rasps, 14mm x 12mm 5mm - 12mm, 7° (+ insert)	1 per height
71-BK-0111	Rasps, 16mm x 11mm 5mm - 12mm, 7° (+ insert)	1 per height

PRECISION SPINE
DAKOTA ACDF™
 STANDALONE CERVICAL SYSTEM


Bone Screws T8

Item No.	Length	Qty/Set
3.5mm Self-Drilling Screws (minor Ø 1.9mm)		
71-DV-3512	12mm	6
71-DV-3514	14mm	6
71-DV-3515	15mm	6
71-DV-3517	17mm	6
4.0mm Self-Drilling Screws (minor Ø 2.4mm)		
71-DV-4012	12mm	6
71-DV-4014	14mm	6
71-DV-4015	15mm	6
71-DV-4017	17mm	6
4.0mm Self-Tapping Screws (minor Ø 2.4mm)		
71-TV-4014	14mm	6
71-TV-4015	15mm	6
71-TV-4017	17mm	6

Plate & Cage Interbody

Item No.	Height*	Qty/Set	Locking Plate Color
14mm x 12mm, 7°			
71-AL-4206	6mm	3	Grey
71-AL-4207	7mm	3	Blue
71-AL-4208	8mm	3	Yellow
71-AL-4209	9mm	2	Green
71-AL-4210	10mm	2	Purple
71-AL-4211	11mm	1	Light Green
71-AL-4212	12mm	1	Pink
16mm x 14mm, 7°			
71-AL-6206	6mm	3	Grey
71-AL-6207	7mm	3	Blue
71-AL-6208	8mm	3	Yellow
71-AL-6209	9mm	2	Green
71-AL-6210	10mm	2	Purple
71-AL-6211	11mm	1	Light Green
71-AL-6212	12mm	1	Pink

*Height does not include total teeth height of 0.6mm



PRECISION SPINE®

www.precisionspineinc.com

