

PRECISION SPINE

SLIMPLICITY[®] PLUS

ACDF SYSTEM

SIMPLE AND VERSATILE

The Slimplicity[®] Plus Anterior Cervical Plate offers one of the slimmest plates available, with an easy-to-use locking mechanism that facilitates visual locking confirmation. Large graft windows have been incorporated to provide an unimpeded graft site and end plate visualization. To optimize bone purchase, a generous 20° screw angulation is provided. The large array of variable and fixed screw options accommodates semi-constrained, constrained and hybrid philosophies.



VARIABLE ANGLED SCREWS

20°



LOW PROFILE DESIGN



ONE STEP LOCKING




PRECISION SPINE[®]

Discover the Difference



SLIMPLICITY® PLUS ANTERIOR CERVICAL PLATING SYSTEM



Low Profile Plate (2mm)
designed to minimize tissue disruption and post-op discomfort

Pre-Contoured
to address patient anatomy



Central Graft Windows
facilitate unimpeded graft visualization

One Step Locking
facilitates secure application

IMPLANT SIZING

Plate Sizing

Lengths*

- 1 Level - 7 sizes, 10-22mm (2mm increments)
- 2 Level - 7 sizes, 25-41mm (3mm increments)
- 3 Level - 8 sizes, 43-64mm (3mm increments)
- 4 Level - 6 sizes, 61-81mm (4mm increments)


*Measured Center Hole to Center Hole

Width 17mm

Waist 14mm

Fixed & Variable Screw Sizes

- 4.0 and 4.5mm Diameters
- 12-18mm Lengths
- (2mm increments)



Single Tray Layout minimizes case complexity


PRECISION SPINE SHURFIT® ACIF 2C HA-CPTI COATED INTERBODY CAGES



- Ti Plasma Coated Surface
- Hydroxyapatite (HA) Coating
- Increased Expulsion Resistance
- Tri-Bond Lock Technology
- Enhanced Imaging Properties

SIZING

- 14mm Wide x 12mm Deep
- 16mm Wide x 14mm Deep
- Heights - 5-12mm
- 5° & 10° Lordotic Options



Precision Spine, Inc.
2050 Executive Drive, Pearl, MS 39208
Customer Service: 1.888.241.4773
Phone: 601.420.4244
Toll Free: 877.780.4370
Fax: 601.420.5501
www.precisionspineinc.com

