


SLIMPLICITY® ANTERIOR CERVICAL PLATING SYSTEM



Low Profile Plate (2mm)
designed to minimize tissue disruption and post-op discomfort

Pre-Contoured
to address patient anatomy

Central Graft Windows
facilitate unimpeded graft visualization

One Step Locking
facilitates secure application

IMPLANT SIZING

Plate Sizing

Lengths*

- 1 Level - 7 sizes, 20-32mm (2mm increments)
- 2 Level - 8 sizes, 37-55mm (3mm increments)
- 3 Level - 8 sizes, 54-75mm (3mm increments)
- 4 level - 6 sizes, 69-89mm (4mm increments)


*Measured end-to-end

Width 17mm

Waist 14mm


Fixed & Variable Screw Sizes

- 4.0 and 4.5mm Diameters
- 12-18mm Lengths
- (2mm increments)



Single Tray Layout minimizes case complexity

PRECISION SPINE SHURFIT® CERVICAL INTERBODY CAGES



- **Large Contact Area**
designed to optimize vertebral body support and minimize risk of subsidence
- **Large Graft Area**
enables excellent biological coverage
- **Trapezoidal Design**
allows for proper anterior placement
- **Tantalum Markers**
strategically placed to aid interbody positioning

SIZING

- 14mm Wide x 12mm Deep
- 16mm Wide x 14mm Deep
- Heights - 5-12mm (1mm increments)
- 5° & 10° Lordotic Options

Precision Spine, Inc.
2050 Executive Drive, Pearl, MS 39208
Customer Service: 1.888.241.4773
Phone: 601.420.4244
Toll Free: 877.780.4370
Fax: 601.420.5501
www.precisionspineinc.com



CE0086