

PRECISION SPINE

# SURELOK™ PC

POSTERIOR CERVICAL SYSTEM

## SIMPLE, VERSATILE SOLUTION


The SureLOK™ PC is a cervicothoracic system that offers a simple and versatile solution to posterior cervicothoracic fixation. The top loading, low profile polyaxial screw design features 80° of angulation to allow for ease of rod insertion with minimal contouring. The system includes easy-to-use Instrumentation, as well as a variety of hook, rod, offset and domino options for the most difficult of cases.



LOW PROFILE



80° ANGULATION



SELF TAPPING, SELF DRILLING




PRECISION SPINE®

Discover the Difference



# SURELOK™ PC POSTERIOR CERVICAL SYSTEM



**Low Profile** Polyaxial Screws provide excellent angulation (80°) in all planes for increased intraoperative flexibility.  
**Self Tapping, Self Drilling** 3.5 & 4.0mm Screws facilitate accurate placement.



**Tulip-to-Tulip** Adjustable Cross-Links provide 3 directions of freedom:

- 30 - 35mm
- 35 - 45mm
- 45 - 65mm

## Full Line of Lateral Offsets, Hooks, Dominoes, Rods and Transition Rods



**Lateral Offsets** available in Open & Closed configurations, 14 & 25mm lengths and a 3.5mm Rod diameter.



**Hooks** available in Straight, Right and Left Offsets, 5 & 6mm sizes.



### Straight and Parallel Domino Sizes

- 3.5 x 3.5mm, 3.5 x 4.5mm
- 3.5 x 5.5mm, 3.5 x 6.25mm

### SureLOK™ PC Implant Sizing

Screws	Cross-Links	Offsets
3.5mm	30 - 35mm	14mm length
4mm	35 - 45mm	25mm length
	45 - 65mm	3.5mm diameter
Hooks	Transition Rods	Rods
5mm (Straight, Right, Left)	3.5 - 5.5mm diameter	3.5mm diameter and 20 - 100mm (10mm increments)
6mm (Straight, Right, Left)	60, 80, 120, 240 & 500mm lengths	120 & 240mm lengths

### Precision Spine, Inc.

2050 Executive Drive, Pearl, MS 39208  
 Customer Service: 1.888.241.4773  
 Phone: 601.420.4244  
 Toll Free: 877.780.4370  
 Fax: 601.420.5501  
[www.precisionspineinc.com](http://www.precisionspineinc.com)

